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Syringe size and selection

- Ensure syringe sizes and models are compatible with the syringe pump (refer to the manufacturer's instructions for use.) Use of incompatible syringe can cause improper pump operation resulting in inaccurate fluid delivery, insufficient occlusion (blockage) sensing, and other potential problems.
- Use the smallest compatible syringe size necessary to deliver the fluid or medication; this is especially important when infusing high risk or life-sustaining medications at low infusion rates (e.g., less than 5 mL per hour, and especially flow rates less than 0.5 mL per hour). Using a larger syringe when infusing at a low rates can lead to inadequate syringe pump performance including delivery inaccuracies, delay of therapy, and delayed generation of occlusion detection alarms. This is due to the increased friction and variable compliance of the syringe plunger tip with larger syringes.

Source : <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm518049.htm>

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2.3.4 Usage considerations

- Syringe pumps have been designed to give optimum performance when placed approximately level with the infusion site. It is not advisable to place the pump well above the infusion site as, even in modern designs, some siphoning can occur in this position.

NB : When the administration set has been connected to the infusion site, the vertical position of the device and giving set in relation to the site should be altered as little as possible. If the pump is raised above the infusion site whilst the liquid is being delivered, it can result in a large bolus being delivered to the patient. (See Section 7 and Appendix 4.)

Source : [UK Medicine and Healthcare Products Regulatory Agency 2013]

Height and Location of the Syringe Pump System

- Ideally, the syringe pump should be level with the distal tip of the catheter (e.g., the site of fluid delivery : if accessing a central line the syringe pump should be at the level of the patient's heart.) If the pump height is raised relative to the distal tip of the catheter (e.g., during patient transport), the increase in height of the syringe pump can result in a temporary increase in fluid delivery or bolus until the flow rate stabilizes. Alternatively, if the pump is lowered relative to the distal tip of the catheter, the decrease in height of the syringe pump may result in a decrease in fluid delivery or under-infusion until the flow rate stabilizes.

Source : <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm518049.htm>

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FREE_FLOW

Although obvious to all, fully opening the roller clamp on a basic gravity infusion set-up will permit all the fluid and any air in the bag to infuse the patient. we seem to abandon this basic logic when using a pump a pump! If the pump is higher than 12 inches above the patient's heart, and the roller clamp is fully open, free-flow is likely occur if :

- The tubing, syringe or cassette aren't carefully located into the housing.
- The pump door is opened or the syringe is removed from it's clamp.
- There is a break/leak in the upper part of the giving pathway or the syringe.

Source : Conti Edu Anesth Crit 2004;4; 81-5, (<http://www.ebme.co.uk/articles/clinical-engineering/46-infusion-devices-training-tutorial>,UKMedicineandHealthcareProductsRegulatoryAgency2013)

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OCCLUSION BOLUS EFFECT

Any blockage/occlusion of the giving pathway causes the downstream line pressure to increase to the pumps occlusion alarm/pressure limit, which causes the tubing to expand/balloon. Clearing the occlusion e.g. opening the tap, roller clamp, kinked tubing will infuse a bolus into the patient. The higher the occlusion alarm/pressure limit is set the larger the bolus.

Source : Conti Edu Anesth Crit 2004;4; 81-5, (<http://www.ebme.co.uk/articles/clinical-engineering/46-infusion-devices-training-tutorial>,UKMedicineandHealthcareProductsRegulatoryAgency2013)